



The MCW 300L chair scale provides an affordable and practical solution to address the challenge of weighing patients with limited mobility or who are unable to stand. MCW combines sturdy construction with functionality, minimising the difficulty of lifting patients and eliminating the need to move heavy health scales. The scale's swivel feet and arm rests add comfort for patients during the weighing process, making it easier for them to get in and out of the chair. Portability is not an issue, as the swivel wheels pivot smoothly, while the scale's built-in battery eliminates the inconvenience of being tethered to an electrical outlet by a power cord. Weight is measured in kilograms and the MCW 300L calculates body mass index (BMI).

## Convenience and accessibility when weighing patients with limited mobility



### Features

- Bright LED display for easy viewing in most environments
- Colour-coded keys facilitate quick recognition of the most frequently used buttons
- Splash-proof housing, keypad and display are easy to clean
- Footrests and armrests pivot for easy patient accessibility
- Wheels swivel, allowing easy movement
- BMI (body mass index) indicator
- Zero-tracking feature ensures display returns to zero reading
- Rechargeable battery included for operation almost anywhere
- AC adapter included

### Applications

- Weighing
- BMI

### About ADAM:

For more than 40 years, professionals worldwide have relied on Adam Equipment for an extensive selection of dependable and affordable weighing equipment. Adam provides the right balance of speed, performance and value.

Model	MCW 300L
<b>Capacity</b>	300kg
<b>Readability</b>	0.05kg
<b>Repeatability (S.D.)</b>	0.1kg
<b>Linearity (+/-)</b>	0.1kg
<b>Weighing Units</b>	kg, lb
<b>Stabilisation Time (sec)</b>	3
<b>Interface</b>	RS-232
<b>Body Mass Index [BMI]</b>	Yes
<b>Height Rod</b>	No
<b>Calibration</b>	External calibration
<b>Display</b>	LED with 22mm -high digits
<b>Power Supply</b>	15VDC 50/60Hz 500mA adapter / Internal rechargeable battery
<b>Operating Temperature</b>	5° to 40°C
<b>Overall Dim.</b>	675x875x950mm (wxdxh)
<b>Net Weight</b>	25kg

## Accessories

Item number	Description
1120011156	ATP thermal printer
3126011263	ATP thermal printer paper
3126011281	ATP thermal printer paper (pack of 10)
3014011014	RS-232 cable M-F
3074010507	RS-232 to USB adapter (must order corresponding RS-232 cable)
600002028	Adam DU data collection program
700660290	Calibration certificate

## Key Accessories



**ATP thermal printer**  
Easy-to-use features combined with quick set-up equals a truly versatile thermal printer. A single push-button feed key and three easy-to-read LED indicators make operation simple.



**Adam DU data collection program**  
Adam DU software connects a balance to the laboratory network, helping to achieve seamless data collection, transfer and analysis. DU features an easy-to-use interface with a main command window, allowing users to quickly collect information from up to eight Adam balances simultaneously using RS-232 or USB connections.

Supplied by